

BARRIER® Neck set

- Effective absorption helping to maintain a dry working area
- Integrated tube holders for instrumentation facilitation
- Generous-sized split sheet for full patient and table coverage



BARRIER® will help you treat more patients to better care

All BARRIER products from Mölnlycke Health Care are designed to offer the optimal combination of protection, functionality and comfort. Our innovative and procedure specific solutions come complete with logistic advantages and expert training and support. These elements together will create a ripple effect of savings in time, effort and costs - helping you treat more patients to better care.

Draping instructions



Place the body drape, labelled 1, in the midline of the patient's chest. Note the direction on the pictogram.



Unfold the body drape towards you, away from you and over the lower part of the body. Remove the release paper.



Fix the adhesive edges to the patient's neck.



Place the head drape, labelled 2, on the patient's head. Note the direction on the pictogram.



Unfold the head drape towards you, away from you and over the patient's head.



Unfold the drape one more time towards the head and remove the release paper.



Fix the adhesive edges.



Draping completed.

Ref. No.	Description	Disp/Trp
697250	Neck set 1 Mayo stand cover, 79x145cm, 65x85cm absorbent 4 Cellulose towels, 18x25cm 2 OP-tapes, 9x49cm 1 Body drape, 175x230cm, adhesive split 20x20cm 1 Head drape, 260x150 cm adhesive split 20x65 cm 1 Table cover, 150x190 cm, 75x190 cm absorbent	7/14

All specified set contents are also available for inclusion in ProcedurePak® trays.

BARRIER surgical drapes

- 100% impermeable materials providing optimal patient safety
- Optimal folding for aseptic application, easy for one person to drape
- User-friendly packaging with dispenser box and peel-open bag for efficient handling

