



BARRIER® Shoulder drape with pouch

- Integrated fluid collection pouch for effective fluid management
- Strong reinforced materials in the critical zone provide reliable quality with high tear resistance
- Soft plastic edges against the wound area give optimal skin adhesion
- 100% impermeable materials providing optimal patient safety
- Optimal folding for aseptic application, easy for one person to drape



Draping instruction



1. Start by placing the stockinette over the hand without unfolding.



2. Note the direction of the pictogram on the first split sheet.



3. Position the split below the arm. Unfold the drape side to side and then towards the patient's feet. Unfold the split section towards the feet.



4. Unfold the stockinette. Fold the absorbent comfort layer of the stockinette in on itself to avoid the risk of microbial transfer. Fix it with the OP-tape.



5. Remove the release paper on the split sheet. Then apply the split as low as possible around the shoulder of the patient.



6. Note the direction of the pictogram on the second split sheet. Unfold the drape side to side and then towards the patient's head. Remove the release papers.



7. Apply the second split overlapping, but not directly on top of, the first split sheet to form a secure barrier.



8. Shape the pouch with the integrated flexi wires. Draping is then completed.

Ref. No.	Description	Disp/Trp
60002	Shoulder set with pouch 1 mayo stand cover 79x145 cm, reinforced 4 cellulose towel 18x25 cm 1 lower shoulder split sheet 230x295 cm, split 10x70 cm 1 stockinette 22x75 cm 2 OP-tape 9x49 cm 1 upper shoulder split sheet with pouch 240x170 cm, split 15x45 cm 1 table cover 150x190 cm, reinforced, 75x190 cm absorbent	4/8

Ref. No.	ProcedurePak® components
70005	Upper shoulder split sheet 240x170 cm, split 15x45 cm
70006	Upper shoulder split sheet with pouch 240x170 cm, split 15x45 cm
70007	Lower shoulder split sheet 230x295 cm, split 10x70 cm

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, SE-402 52 Göteborg, Sweden.
 T +46 31 722 30 00. F +46 31 722 34 00. www.molnlycke.com

The Mölnlycke Health Care, BARRIER®, and ProcedurePak® names and respective logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies. Copyright (2012)