

Extra Large Leggings

- Integrated leggings that provide full coverage for all lithotomy equipment
- Built in adhesive allows user to adjust legging to their own needs
- 100% impermeable
- Fast draping, easy for one person to use
- Compatible with “banana” boots and stirrups



You know how you must maintain the sterile field when operating on patients in the Lithotomy Position? And how sometimes the size of the foot holder can make this difficult? What we have done is to increase the size of our leggings so that this is no longer a problem. We have also provided two adhesive areas on each legging so that you can safely secure any excess drape material out of the way if needed. This means that, however your patient's leg is supported, the new legging will do the job by providing complete and secure coverage.



01 This is an example of a banana boot, illustrating the need.



02 Banana boot with a patient in the lithotomy position.



03 Draping.



04 Removing backing paper.



05 Reducing drape height.



06 Reducing drape length.

Sterile Pack Code	BNS Code of main drape	Description	Previous design		New design	
			Leg length (cm)	Leg width out/in (cm)	Leg length (cm)	Leg width out/in (cm)
888222-99	965520	TUR Set	120	87	135	100
888225-99	965525	TUR Set	120	87	135	100
888224-99	965522	TUR Set	120	87	135	100
696600-01	956966	Gynaecology/Laparoscopy Set	125	25/75	135	100/100
888112-66	965523	Gynaecology Set	125	100	135	100
888113-99	965523	Gynaecology Set	125	100	135	100
888212-99	965523	Gynaecology Set	125	100	135	100
694240-99	944200	Laparoscopy Set Abdo-Perineal	120	87	135	100
694241-00	944200	Laparoscopy Set Abdo-Perineal	120	87	135	100
694245-99	944205	Laparoscopy Set	120	87	135	100
696810-70	969681	Gynaecology Set	120	75/100	135	100/100

