

The Gynaecology Range


- 100 % impermeable all over
- Fast draping, easy for one person to drape
- Absorbent working area
- Self-adhesive edges around the vaginal area
- Fluid-control through integrated collection pouch



The Gynaecology Range


Gynaecology Set / Ref. No. 696700	Size (cm)
4 cellulose towel	19 x 25
2 op-towel	75 x 90
1 adhesive towel	50 x 50
2 op-tape	9 x 49
1 lithotomy drape with leggings	150 x 200
1 table cover	150 x 190

Pack Count: 7/14




Gynaecology Set Ref. No. 888112 / 888212	Size (cm)
1 mayo stand cover (888112 only)	79 x 145
2 cellulose towel	19 x 25
1 tube holder	9 x 18
1 gyn drape with pouch	60 x 116
1 base drape with leggings	200/270 x 215
1 table cover	150 x 190

Pack Count: 7/14




Gynaecology Set / Ref. No 888113	Size (cm)
1 mayo stand cover	79 x 145
2 cellulose towel	19 x 25
1 tube holder	9 x 18
1 gyn drape	60 x 116
1 base drape with leggings	200/270 x 215
1 table cover	150 x 190

Pack Count: 7/14



Uro/Gynaecology Set Ref. No. 905020	Size (cm)
1 mayo stand cover	79 x 145
1 under buttocks drape	75 x 100
4 cellulose towel	19 x 25
1 uro/gynaecology set	200/355 x 295
1 table cover	150 x 190


Pack Count: 5/10



The Gynaecology Range


Gynaecology Set / Ref. No. 925984	Size (cm)
1 tube holder velcro	2,5 x 30
1 mayo stand cover	79 x 145
4 cellulose towel	19 x 25
1 under buttocks drape	75 x 100
1 gyn drape	200/270 x 280
1 table cover	150 x 190

Pack Count: 5/10




Gynaecology / Cystoscopy Set Ref. No. 696110	Size (cm)
2 leggings	75 x 120
1 aperture drape	75 x 175
1 table cover	150 x 190

Pack Count: 10/20




Gynaecology / Cystoscopy Set Ref. No. 902496	Size (cm)
1 utility drape	75 x 80
2 leggings	75 x 120
1 cystoscopy drape	75 x 175
1 table cover	150 x 190

Pack Count: 9/18



Gynaecology / Laparoscopy Set Ref. No. 696600	Size (cm)
2 cellulose towel	19 x 25
1 op-tape	9 x 49
1 adhesive towel	75 x 75
1 laparoscopy drape, with leggings	250/270 x 310
1 table cover	150 x 190

Pack Count: 7/14





Mölnlycke Health Care AB, Box 13080, SE-402 52 Göteborg, Sweden. T +46 31 722 30 00. F +46 31 722 34 00. www.molnlycke.com
The Mölnlycke Health Care name and logo and BARRIER® are registered trademarks of Mölnlycke Health Care AB.